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The screenshot displays the Recipe Studio interface for a recipe named "AMOX_WITH_TRAINS". The main workspace shows a flowchart with the following steps:

- BEGIN PROCEDURE EXECUTION
- Parallel steps: UP_FILL_MIT:1, UP_FILL_MIT:2, UP_FILL_MIT:3
- WAIT UNTIL INGREDIENT TANKS FILLED
- Parallel steps: UP_RDY_TO_XFR:1, UP_CHG_MPT_PXFR_FRM:1
- PREP TANK SELECTED AND READY TO START TRANSFER

Key interface elements are highlighted with yellow boxes:

- 工具栏 (Toolbars):** Located at the top, containing menus like Home, Diagram, View, and various icons for editing and execution.
- 配方步骤 (Recipe Steps):** A yellow box highlighting the central flowchart area.
- 参数定义 (Parameter Definition):** A yellow box highlighting the parameter table at the bottom.

The parameter table is as follows:

Parameter	AMOX_HIGH_CSL	AMOX_HIGH_CSL_LAC	AMOX_HIGH_LAC	AMOX_STANDARD
AGITATE_TIME	17	14	14	<10>
CSL_AMOUNT	144.000000	144.000000	50.000000	50.000000
LAC_AMOUNT	20.000000	100.000000	100.000000	50.000000
MIT1_LEVEL	<700.000000>	<700.000000>	450.000000	250.000000
MIT2_LEVEL	450.000000	<700.000000>	<700.000000>	250.000000
MIT3_LEVEL	450.000000	450.000000	450.000000	250.000000
RSD_AMOUNT	20.000000	18.000000	14.000000	20.000000

On the left, a tree view shows the recipe structure, and a table below it lists parameters with their values, types, and categories.

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Select Recipe

Unit selections for user: ADMINISTRATOR Change User

Error:

General Info

Batch ID: DEMO_UNIT_ALIAS

Scale: Min: 0.00, Value: 100.00, Max: 100.00, EU: %, Scale: 100.00 %

Recipe

Name: AMOX_WITH_TRAINS

Version: 2 Details

Description: AMOXICILLIN FERMENTATION MEDIA PREPARATION WITH C

Equipment Train

Name: OP_SELECT_AREA1

Description: Operator Selects Units for Area 1

Unit Aliasing

Alias Name	Unit Name	Selection Mode
#NGRED_TANK1#	{Defer}	DEFER
#NGRED_TANK2#	{Defer}	DEFER
#NGRED_TANK3#	{Defer}	DEFER
#WFI_TANK#	WFI-1400	DEFER
#PREP_TANK#		FROM_POLICY

Use upstream/downstream filtering of unit modules

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Select Recipe

Unit selections for user: ADMINISTRATOR

Error:

General Info

Batch ID: DEMO_BATCH

Scale

Min	Value	Max	EU	Scale
0.00	100.00	100.00	%	100.00 %

Recipe

Name: AMOX_WITH_TRAINS

Version: 2

Description: AMOXICILLIN FERMENTATION MEDIA PREPARATION WITH C

Equipment Train

Name: 1500_TRAIN

Description: No Operator Selection, 1500 Prep Tank

Unit Aliasing

Alias Name	Unit Name	Selection Mode
#NGRED_TANK	MIT-1100	DEFER
#NGRED_TANK	MIT-1200	DEFER
#NGRED_TANK	MIT-1300	DEFER
#WFI_TANK#	WFIT-1400	DEFER
#PREP_TANK#		FROM_POLICY

Use upstream/downstream filtering of unit modules

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
Batch ID	Description	Value
AMOX_BATCH_123	Ready to begin making media solution?	
AMOX_BATCH_456	Select Equipment For UP_CHG_MPT_P1/P2_FR	

Batch ID	Comment	Operator
AMOX_BATCH_456	This batch running on train 3 due to maint	Operator:ADMINISTR

Batch Event Journal						
2/22/2009 10:01:11 PM	AMOX_BATCH_456	A1_AMOX		Header		
2/22/2009 10:01:42 PM	AMOX_BATCH_456	A1_AMOX_MAN_MEDIA_PREP	Comment			This batch running on train 3 due to maintenance.
	AMOX_BATCH_456	A1_AMOX_MAN_MEDIA_PREP	Step Activated	Step Activity		Initial Step
	AMOX_BATCH_456	A1_AMOX_MAN_MEDIA_PREP	Procedure Started	System Message		0

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DeltaV Operate (Run)

Module: Main: Bio_Ovw_1000 Username: ADMINISTRATOR 3:23:01 PM

Area 1000 Overview

Request Dialog

Do you want to hold 'AMOX_BATCH_456'?

Comment: Hold batch while performing filter maintenance.

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


- “Rebind Units” gives operators the agility to switch units right down to the last minute – even if the batch is already running
- Minimize downtime and keep batches running when unforeseen circumstances arise in the plant
- Easy to correct invalid operator equipment selections

Increased operational flexibility!

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