





## 配液系统基本要求



- Water supply
- Material dosing
- Vent/Vacuum/Sparge
- Mix
- Heat/Cooling
- Product transfer
- Filtration
- Sample
- Storage

5







WINA TECH


## 多路冷点回换热器的原因及优点


原因：


由于冷点比较多，如果每次用水结束后直接回主管道的话，会造成主管道的温度下降，这点在GMP的角度上会有很大问题，因为大面积的冷水会造成很大风险


优点：

把这些冷水直接回到末端换热器的前端，这样的话 这些冷水直接通过换热器进行加热，从而使风险降到最低









用点连接3



**WINA TECH**

**Loading cell/Balance**

Accuracy 0.1%  
Economy  
Replace level meter

Hose connection  
Agitator running impact  
Installation difficulty  
Calibration online



12



- ❖ Manual valve
- ❖ Control valve
- ❖ On-off auto valve
- ❖ On-off auto valve with small by-pass valve




**Material dosing**

**WINA TECH**

**Manual dosing into man hole**

- Hopper gravity
- Vacuum dosing system
- Peristaltic pump
- Powder dispersion

14



WINA TECH

Powder dispersion

Recirculation

16



**Cooling/Heating**


WINA TECH

Heating and cooling by jacket

Cooling by jacket

TCM

17





**Product transfer**

WINA TECH

Sterile PA/N2 transfer

Lobe pump

Peristaltic pump

Centrifugal pump

Quattroflow diaphragm pump

19




**Single filter**


**Redundant filter**

**Liquid filter integrity test in-line after SIP and product**

**Filter layout**

21








The image shows a close-up, slightly blurred view of various pieces of laboratory glassware, such as test tubes and flasks, submerged in a blue liquid. In the top right corner of the collage, there is a logo for "WINA TECH" featuring a stylized blue wave graphic above the word "WINA TECH". Below the collage, there is a white rectangular area containing text. At the very bottom of the slide, there is a thin horizontal bar with a blue gradient.


**Sterile Sample**


Mixing temperature control


N2 protection

24









**WINA TECH**

## Spray Ball

| Type            | Advantages   | Disadvantages  | Illustration |
|-----------------|--|--|--------------|
| Low pressure    | <ul style="list-style-type: none"> <li>Low cost</li> <li>Wide range available</li> <li>Rapid heat input</li> <li>No moving parts</li> <li>Low maintenance</li> <li>Copes well with internal obstructions</li> <li>Flowrate integrates well with line cleaning</li> <li>Easy to confirm correct operation</li> </ul>  | <ul style="list-style-type: none"> <li>High flowrate</li> <li>High effluent and water load</li> <li>No mechanical effect</li> <li>Relies entirely on chemical cleaning</li> <li>Limited spray diameter</li> <li>High pressure vapourise sprays</li> <li>Need multiple spray heads for larger vessels</li> <li>Need larger CIP pipework and scavange pumps</li> </ul> |              |
| High pressure   | <ul style="list-style-type: none"> <li>Low flowrate</li> <li>Lowest effluent load</li> <li>Removes heavy soil by high-pressure jets</li> <li>Lowest CIP pipework diameters</li> <li>Smaller scavange pumps</li> <li>Excellent spray diameter range</li> <li>High-pressure mechanical cleaning effect</li> <li>Single head can clean very large vessel</li> </ul> | <ul style="list-style-type: none"> <li>High cost</li> <li>Lower heat input rate</li> <li>Moving parts</li> <li>Higher maintenance</li> <li>Can jam</li> <li>Needs flow and motion instrumentation to automatically confirm operation</li> <li>Internal obstructions break up jets</li> <li>Longer cleaning sequence</li> </ul>                                       |              |
| Medium pressure | <ul style="list-style-type: none"> <li>Medium cost</li> <li>Wide range available</li> <li>Moderate heat input</li> <li>Medium flowrate</li> <li>Lower effluent load</li> <li>Lower CIP pipework diameters</li> <li>Smaller scavange pumps</li> <li>Good spray diameter range</li> <li>Some mechanical cleaning effect</li> </ul>                                 | <ul style="list-style-type: none"> <li>Need multiple spray heads for larger vessels</li> <li>Moving parts</li> <li>Higher maintenance, and can jam</li> <li>Need flow and motion instrumentation to automatically confirm operation</li> <li>Internal obstructions break up jets</li> </ul>  |              |

**TANKO® RT**由旋轉洗球 TANKO® RPB-E (A), 液壓管 (B) 及氣動驅動裝置 (C) 組成一體


29



Bio-tech CIP Step:

- ① Drain
- ② Leakage test
- ③ PW Pre-washing (cold)
- ④ Detergent(1~3% concentration) heating
- ⑤ Detergent washing and recirculation
- ⑥ Air blow to reclaim detergent
- ⑦ PW rinse
- ⑧ WFI final rinse
- ⑨ Air blow to dry

CIP flow rate and pressure variation control within 15%



**FIGURE 9** Filter housings absent nozzles in the dome may be cleaned with a simple disc distributor. This device may be left in place for all use of the filter and will distribute flush, wash, and rinse solutions to the housing dome and sidewall at any range of flow and pressure required to meet other CIP design criteria.




## Step:

- ①Drain
- ②Leakage test
- ③Vacuum
- ④Clean steam heat up
- ⑤Sterilization
- ⑥Air blow cooling down
- ⑦Cold medium cooling down
- ⑧Pressurized



## Typical SIP drip tube installation detail





WINA TECH

## Equipment Selection

## 设备选型



### 罐底阀 Tank Bottom Valve

- 材质: 316L
- 标准: ASME BPE, 特殊阀体整块钢加工成, EP Ra<0.4um
- 膜片: 双层膜片EPDM/PTFE
- 优点: 焊在容器的最底部, 从而可以达到最佳的排空、清洗、消毒的效果。阀的密封堰已经尽可能靠近罐底壁, 从而达到零死角效果。


阀体由整钢加工而成（无焊接的结构）





### 磁力搅拌器 Magnetic Mixer

- 接触物料材质: 316L
- 表面处理: EP Ra<0.4um
- 其他: 碳化硅轴侧, 不能干转。  
设计转速50-400rpm, 叶轮锻造,  
同轴异步电机。含测数器, 防护  
等级IP55





### 罐体视镜 Sight Glass

- 型号: MV80+SWI+USL01
- 材质: SS316L
- 观察孔: soda lime glass
- 操作压力: 6bar
- 温度: 最高280°C



MV/ME系列，  
可满足光源和  
观察口的用途



MV+USL01系  
列

## Equipment Spec 设备选型



### 回流泵 pump

- 物料接触部分: SS316L
- 抛光: MP Ra 0.6um
- 密封: 单端机械密封, SiC/SiC/EPDM

耐温130°C



## Equipment Spec 设备选型



### 管壳式换热器

- 密封圈：采用符合 FDA 标准的 O 形密封圈。
- 设计形式：管层设计为双层壳层走蒸汽 物料走管层






## Equipment Spec 设备选型


### 疏水器

型号: BTM7洁净型热静力疏水器,  
卡箍连接  
安装位置: CIP灭菌管路


BTM7  
卫生卡箍式接头



型号: MST21H热静力疏水器, 丝扣  
连接  
安装位置: SIP排放



型号: FT14浮球式疏水器, 法兰  
连接  
安装位置: PS冷凝水排水管路





## 仪表Instrument

### CIP系统专用电导率仪表 Conductivity meter

- 测量范围 Range of testing: 600-2000us/cm
- 输出 Output: 0-20mA
- 精度 Precision : ±1.0%
- 温度上限 The temperature upper limit : 155°C
- 连接 Connection: Clamp connection 卫生型卡箍连接
- 测定点 Testing: 在化学添加剂混合后打出的管道上




## Equipment Spec 设备选型



### 流量计 Flow meter

- 质量流量计
- 材质：316L不锈钢
- 测试范围：0-45m<sup>3</sup>/h
- 介质温度：100°C,
- 输出：4-20mA
- 精度：0.5%
- 表面抛光：MP Ra<0.4μm
- 连接方式：卫生型卡箍连接
- 安装位置：收获系统







## System Function Test

- ✓ System Software Edition Check 检查系统软件版本
- ✓ Password Test 密码测试
- ✓ System Auto/Man/stop run 系统自动 手动及停止
- ✓ CIP pump running CIP 泵运行
- ✓ Tank level 罐子液位
- ✓ CIP/SIP 测试
- ✓ Filter Integrity Test 过滤芯完整性检测
- ✓ Cooling down test 冷却测试 (City Water)
- ✓ Vent Filter Heat Jacket Function Test 呼吸器夹套加热功能测试
- ✓ System emergence stop test 系统紧急停止测试
- ✓ Agitation Test(VFD) 搅拌器调速测试
- ✓ Load Cell Test(Raw material adding test)
- ✓ Flow meter Test(Transfer Test)




**Test Facilities 测试设施**


**WINA TECH**

Test Facilities Include:

测试设施包括:

- RO Water  
RO水
- Plant Steam  
工厂蒸汽
- Clean Steam  
洁净蒸汽
- Clean Dry Air  
清洁干燥空气
- Vessel Test Bay  
容器试验台
- Sprayball testing  
喷淋测试
- Calibrated inspection equipment includes:  
Micrometers, Profilometers, Ferritescope,  
Digital Protractors, Positive Material  
Identification (PMI).  
校准检测设备包括: 千分尺, 表面粗糙度仪, 铁离子检测仪, 数显角度仪, 材料检测仪 (PMI)。







**SAT**  
现场验收测试

**WINA TECH**

## Function Test系统功能测试

- 密码测试
- 泄漏测试
- CIP测试
- 吹扫测试
- 灭菌测试
- 过滤器湿润
- 完整性测试
- 称重校验
- 配制测试
- PH调节
- 搅拌测试
- 定容测试
- 物料传输测试
- 报警测试
- Etc.

**cGMP Perspectives**

**WINA TECH**


## cGMP Perspectives

**For single line**

- Material
- Deadleg
- Drainability
- Welding
- Calibration
- Marking & Labeling
- Insulation
- Pickling & Passivation
- Decontamination
- Valve Angles
- Air Break
- Filter Integrity
- Etc.

**For System**

- Surface Finish test
- Sprayball Coverage Test
- GAMP5 & Part 11
- Minimum Flow Rate
- Wiring Test
- System Alarm
- Interlock
- P&ID Verification
- Critical Process Parameter
- Critical Quality Attributes
- Etc.





**Material Inspection**

Qualified Materials

**Sample Installation**

Material Certificate

Surface Finish Report

Random Check

3<sup>rd</sup> Party Check

This slide is titled "Material Inspection". It features a purple header bar with the text "Qualified Materials". Below this is a light blue rectangular area containing the heading "**Sample Installation**". At the bottom of this area are four small square icons with labels: "Material Certificate" (blue), "Surface Finish Report" (purple), "Random Check" (dark blue), and "3<sup>rd</sup> Party Check" (teal). The background of the slide has abstract blue and yellow wavy patterns. In the top right corner, there is a logo for "WINA TECH" featuring a stylized blue wave icon above the company name.



**WINA TECH**


### Equipment & Welder certificate 施工设备及焊工资质




**Imported endoscope**  
进口手持式内窥镜3pcs



**spectrum analyzer**  
手持式光谱分析仪1pcs



**AWS certificate**  
AWS 证书



**SGS certificated welder**  
国际第三方认证焊工

The slide features a blue-toned background image of a welding torch and sparks. In the top right corner is the WINA TECH logo, which includes a stylized blue wave icon above the text "WINA TECH". The main title "Welding Inspection" is centered in a large, bold, white font. Below the title is a purple rounded rectangle containing the text "Qualified Welding". A large light blue rectangular area contains the text "3<sup>rd</sup> Party Inspection". At the bottom, there are four smaller blue rectangular boxes with white text: "WPS & Welder Qualification", "Purging Gas Quality", "Boroescope", and "Welding Log & Map".

**Qualified Welding**

**3<sup>rd</sup> Party Inspection**

WPS & Welder Qualification

Purging Gas Quality

Boroescope

Welding Log & Map





The slide features a blue and yellow abstract background at the top. In the top right corner is a logo for "WINA TECH" with a stylized blue wave graphic above the text. Below the logo, the text "Welding Quality Control 焊接质量控制" is displayed in white. A large blue downward-pointing arrow is positioned in the center. To its left is a blue rounded rectangle containing the text "文件记录". To the right of the arrow are five white rectangular boxes, each containing a piece of welding quality control information: "焊点图及焊点编号", "每日焊样记录", "每日焊接记录", "内窥镜记录", and "第三方抽检记录". The bottom of the slide has a solid blue horizontal bar.

WINA TECH

Welding Quality Control 焊接质量控制

文件记录

焊点图及焊点编号

每日焊样记录

每日焊接记录

内窥镜记录

第三方抽检记录

**WINA TECH**

**Welding Quality Control**  
焊接质量控制

预制配管及准备

The Approved Isometric Drawing  
被批准的轴测图

Pipe cutting and facing  
管道切割及切口平口

High purity Ar and PTFE pipe  
合资高纯氩气及PTFE焊接气管

Spot Welding  
焊接管对口及三  
点式点焊

The diagram shows a vertical sequence of four main steps in blue boxes, each with a corresponding image:

- Pre-fabricated piping preparation:** Shows a complex blue and yellow isometric piping layout.
- The Approved Isometric Drawing:** Shows a red and black pipe cutting and facing machine.
- High purity Ar and PTFE pipe:** Shows a large stainless steel gas cylinder and a roll of white PTFE tape.
- Spot Welding:** Shows a worker wearing a yellow hard hat and safety glasses performing a weld on a metal pipe.

**WINA TECH**

**Welding Quality Control**  
焊接质量控制

The diagram illustrates the setup of an oxygen content test machine. It shows an Argon Dewar (99.999%) connected via a tube to a Regulator, which then connects to a Flowmeter and a Filter (0.02 Micron). The filtered gas then flows through a tube to a Weld, where it is used for Purge Gas. An Oxygen Analyzer is connected to the system. The chart on the right shows eight photographs (A-H) of welds with increasing oxidation levels corresponding to oxygen concentrations of 21 ppm, 31 ppm, 100 ppm, 200 ppm, 1,000 ppm, and 10,000 ppm.

| Concentration (ppm) | Color Description |
|---------------------|-------------------|
| 21 ppm              | A                 |
| 31 ppm              | B                 |
| 100 ppm             | C                 |
| 200 ppm             | D                 |
| 1,000 ppm           | E                 |
| 10,000 ppm          | F                 |
| - ppm               | G                 |
| - ppm               | H                 |

**Oxygen content test machine**  
含氧量测试设备

**Weld Oxidation Color**  
焊缝氧化颜色

Boroscope ratio在焊缝检查中，内窥镜检查比例：  
Module模块部分 检查比例为100%  
Site现场部分 Pre-fabricate 100%预制检查比例为100%，site piping 30-50%现场管道安装  
检查比例为30~50%

**Pickling & Passivation 酸洗钝化**

**Standard Material 标准材料:**

波洱Deconex CIP Power-x 浓缩高效碱性清洗液  
波洱Deconex CIP Protect 浓缩液体钝化清洗液

**Standard Procedure 标准清洗流程:**

Refer to the next page 请见下页

**CERTIFICATE OF GMP COMPLIANCE**

We certify herewith that the company **Borer Chemie AG** with its site of **Gewerbestrasse 13, 4528 Zuchwil, Switzerland**, has been duly authorized to manufacture and distribute medicinal products, the manufacturing licence excluding sterile products and including following dosage forms:  
 - Liquid dosage forms (restricted to disinfectants)  
 that the finished medicinal products put on the market in Switzerland by the company are subject to approval and authorisation by our agency;  
 that the company is keeping the required level for good practices in the manufacture of pharmaceutical products and active pharmaceutical ingredients according to the Swiss regulations in force. These regulations are in accordance with the requirements for good practices in the manufacture and quality control of the Pharmaceutical Inspection Convention /Co-operation Scheme (PIC/S) and the Directives of the European Commission;  
 that the manufacturing plant of the company is subject to official periodic inspections; the last regular inspection was conducted November 16-17, 2010;  
 that the requirements regarding manufacture and quality control for pharmaceutical products and active pharmaceutical ingredients for export are identical to those applicable to products sold in Switzerland.

Berne, Januar 12, 2011  
No. 11-051

Swissmedic, Swiss Agency for Therapeutic Products  
Dr. Georges Meléguer


Swissmedic, Swiss Agency for Therapeutic Products  
Gewerbestrasse 13 • Postfach • CH-8006 Zürich • Tel.: +41 31 522 52 00 • Fax: +41 31 522 52 10



**Pickling & Passivation Standard Procedure**  
酸洗钝化标准清洗流程

**WINA TECH**


| No.<br>序号 | Steps<br>操作步骤   | Media<br>介质参数   | Temperature<br>温度              | Retention time<br>保持时间                                 |
|-----------|---|---|--------------------------------|--|
| 1         | Pre-wash 预洗   | Use Drinking water to rinse<br>饮用水冲洗                    | At room temperature 室温         | 1-2min   |
| 2         | Step one 清洗第一步:<br>To remove the organism<br>去除油脂等有机物 | 2%, CIP Power X solution 溶液                             | 60-70 °C                       | Re-circulation 循环 15-30min                             |
| 3         | Step one<br>漂洗第一步                                     | Use PW to rinse<br>纯化水冲洗                                | At room temperature 室温         | 1-2min   |
| 4         | Step two<br>漂洗第二步                                     | Use PW to rinse<br>纯化水冲洗                                | At room temperature 室温         | 1-2min   |
| 5         | Step two : Passivation<br>漂洗第二步：钝化                    | 3-5%, CIP Protect solution 溶液                           | 65-75 °C                       | Re-circulation 循环 60min                                |
| 6         | Step three:<br>漂洗第三步                                  | Use PW to rinse<br>纯化水冲洗                                | At room temperature 室温         | 1-2min   |
| 7         | Step four:<br>漂洗第四步                                   | Use PW to rinse<br>纯化水冲洗                                | At room temperature 室温         | 1-2min   |
| 8         | Step five:<br>漂洗第五步                                   | Use PW to rinse<br>纯化水冲洗                                | At room temperature 室温         | 1-2min   |
| 9         | Sanitization<br>消毒或者灭菌                                | Medical sanitization or PS sanitization<br>化学消毒或者高温蒸汽灭菌 | Room tem. or 室温或者<br>125-140°C | According to the requirement<br>按照需要 (不建议使用含氯类<br>消毒剂) |






❖ E 阀门安装角度检查 Installation valve and angle inspection

在水平管路上安装的隔膜阀，生产商必须明确需要安装的角度（如果不能以其他方式排水）。测量的误差不需考虑。





F 空气隔断的检查 Air break inspection





# Questions? Thank You !

More information please contact  
William Wang  
[wiliamw@winatech.com.cn](mailto:wiliamw@winatech.com.cn)