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| Concentration [at %] Sampling location | Cr / Fe | Cr ox. / Fe ox. |
|--|---------|-----------------|
| Surface condition: brushed without passivation | 0.3 | 0.3 |
| Surface condition: brushed HNO ₃ passivation | 0.8 | 1.0 |
| Surface condition: brushed Electropolishing | 0.8 | 1.0 |
| Surface condition: electropolished HNO ₃ passivation | 1.2 | 1.9 |

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Substrates stainless steel 1.4404

| Finish of surface | electropolished | | ground | | handpolished | |
|------------------------------------|-----------------|---------------|--------|---------------|--------------|---------------|
| concentration [at %] sample | Cr/Fe | Cr ox./Fe ox. | Cr/Fe | Cr ox./Fe ox. | Cr/Fe | Cr ox./Fe ox. |
| process 1 original process | 2.2 | 4.3 | 6.0 | 7.4 | 3.7 | 5.3 |
| process 2 shorter contact times | 1.8 | 2.8 | 7.0 | 7.0 | 2.5 | 3.3 |
| reference value | 1.1 | 1.6 | 0.4 | 0.4 | 0.3 | 0.3 |

| | 1.4301 X5CrNi1810 brushed surfaces | | 1.4034 X46Cr13 brushed surfaces | | 1.4021 X20Cr13 brushed surfaces | |
|------------------------------------|--|-----------------|---------------------------------------|-----------------|---------------------------------------|-----------------|
| concentration [at %] sample | Cr / Fe | Cr ox. / Fe ox. | Cr / Fe | Cr ox. / Fe ox. | Cr / Fe | Cr ox. / Fe ox. |
| process 1 original process | 2.0 | 3.8 | 2.2 | 3.0 | 1.6 | 2.4 |
| process 2 shorter contact times | 1.7 | 2.8 | 3.0 | 3.6 | 2.0 | 2.5 |
| reference | 0.3 | 0.3 | 0.2 | 0.2 | 0.1 | 0.1 |

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**deconex®
DEROUGE /
DEROUGE P**

Derouging system at neutral pH value

Rapidly and effectively removes rouge deposits from stainless steel water systems and process vessels

Application

The deconex DEROUGE system rapidly and effectively removes rouge from the surfaces of AISI 316L, AISI 316 TI and AISI 904L stainless steel. The new process works in neutral pH conditions and therefore meets the ever more stringent requirements for plant maintenance and environmental protection.

The deconex DEROUGE system is a controlled process and is therefore suitable for derouging production equipment in the following manufacturing sectors:

- Pharmaceuticals
- Biotechnology
- Cosmetics
- Purified water production for medical uses

The deconex DEROUGE system is especially suited to the derouging of:

- Water for injection (WFI) systems (loops and storage tanks)
- Clean steam systems
- Process vessels (e.g. in vaccine production)
- Autoclaves

Properties

As the natural passivated layer on stainless steel surfaces ages it can form a rust-red deposit known as rouging. Rouging is generally considered to be unacceptable because there is the risk that particles of the deposit will detach and enter circulation.

Rouge deposits also make it difficult to control the process parameters of systems where a specified surface state must be maintained.

Your solution for removing rouge: the deconex DEROUGE system (patent applied for) is an innovative, neutral pH derouging method which achieves excellent and rapid results.

The system consists of deconex DEROUGE liquid and deconex DEROUGE P powder. These two components are mixed shortly before use.

The advantages of using the deconex DEROUGE systems are:

| Advantages | Benefits to you |
|----------------------------------|---|
| Short process times | Shorter downtimes, longer duty cycles for production equipment |
| Neutral pH process | Solutions and rinsing water do not require neutralising. Waste disposal is simple and cost-effective. |
| Safer to use | High concentrations of aggressive acids are not required. |
| Simple in-process control | Monitoring of derouging solution effectiveness does not require expensive analytical equipment. |
| Rapid results | Short derouging cycle times. Improved efficiency. |
| Tried-and-tested process | Safer to use thanks to our process descriptions. This permits easy documenting of derouging operations. |
| Highly effective | The lower concentrations required mean that it is more cost-effective to use. |
| Excellent material compatibility | Longer life cycles for process plant and equipment. |
| Does not contain halogenides | No chemical attack of surfaces. |
| Validated analytical methods | Allows a residue-free process. A well established analytical procedure enables rapid transfer of the derouging process. |

DEROUGE P

Complexing agent and corrosion inhibitor
Reducing agent

Iron < 0.5 %
Depends on type and characteristics

1,4421, 1,4404, 1,4426, 1,4571, 1,4539
UDF, EPO, silicone, PEEK

pH value
Density
Appearance
pH value
Appearance

1

Neutralization

A neutralization is not necessary. Process and rinsing water can be disposed of without the need for further treatment. Dispose of process and rinsing water in accordance with your local waste water regulations.

Residue analysis/cleaning validation


Analysis methods for testing for the presence of process chemical residues are available.

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
Treatment with new deconex DEROUGE system. Before and after 1 hour.

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

| | Strong acid | Soft acid | pH-neutral |
|-----------------------------|------------------|-----------|-------------------|
| Efficacy | - | - | ++ |
| Health, Safety, Environment | -- | - | + |
| Time | - | - | ++ |
| Complexity of Operation | +/- | +/- | +/- |
| Cost | | | |
| Material compatibility | +/- | + | ++ |
| Application example | 30 %, 70 °C, 4 h | | 1.3 %, 70 °C, 4 h |

Before and After
Derouging on a piece of piping

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
Sample image (256 μm \times 256 μm) before (left) and after (right) derouging; right: confocal image

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
Before/after pH-neutral derouging treatment

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Corrected for surface waviness error

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