# IHOLLAND®

T A B L E T T I N G S C I E N C E

**IHOLLAND**  
Meadow Lane, Long Eaton  
Nottingham, NG10 2ED  
England  
Tel: +44(0)115 952 6153  
Fax: +44(0)115 952 1188  
Email: [info@iholland.co.uk](mailto:info@iholland.co.uk)  
Web: [www.iholland.co.uk](http://www.iholland.co.uk)

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**Maintain Good Punch Tip Face Finish**  
I Holland PharmaCare<sup>®</sup>


Following a maintenance procedure can ensure that all essential steps are carried out in order to prolong tooling life.

Assessing and polishing your tooling are two crucial steps to ensure that any defects found on the tooling can be dealt with at early stages.

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
**Steel Selection**  
I Holland PharmaGrade®

I Holland's range of steels achieves an optimum balance of a variety of characteristics. One being strength and its ability to withstand an applied stress without failure.


**Suitable Tablet Profile**  
I Holland Design

Various tablet shapes are unable to withstand a high compaction force due to the weight distribution, e.g. deep cup profiles will have a softer core whilst flat bevelled edges have a more even distribution of compaction force.

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
Suitable Dwell Time  
XDF™ (eXtended Dwell Flat)

I Holland's innovative patented elliptical head form can increase dwell time by up to 50%. If the sticking issue is due to the amount of time the granule gets compacted, XDF is the solution.

XDF allows for more dwell than a D-type punch on a B-type tool, allowing your customers to solve dwell problems without upsizing their punches or making a significant investment on a new press.

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I Holland (Shanghai) Limited